



Anbessa Imports

ANBESSA MEDICAL · STRATEGIC SOURCING & SUPPLIES

12-Lead ECG Machine

CARDIAC CARE



Interpretive 12-lead electrocardiograph with auto-analysis, built-in thermal printer, and USB/WiFi connectivity.

Price: TBC — request a quote

Key Benefits

- Automatic interpretation to support rapid clinical decisions
- Built-in thermal printer and 1,000+ record storage
- USB and WiFi transfer for record keeping

Specifications

Channels	12-lead simultaneous
Display	7" color LCD
Sensitivity	2.5–40 mm/mV
Print Speed	25/50 mm/s
Storage	1000+ records

About Anbessa Imports

Anbessa Imports is the sister company of Anbessa Travel, an Ethiopian business known for its sustainability certification, gold level taxpayer status and 13 years of exceptional customer service. We replicate that same heritage in medical equipment: responsibly sourced from vetted manufacturers, ethically managed through every stage of procurement and compliance, and delivered with installation, training and after-sales support across Ethiopia.

- Sister company of Anbessa Travel
- Sustainability certified, gold level taxpayer
- 13 years of exceptional customer service
- Responsibly sourced and ethically managed

Anbessa Imports — Anbessa Medical

Bole Sub-City, behind the Japanese Embassy, on the road between Bole Sub-City HQ and Egzabherab Church, in front of Chiko Restaurant, Addis Ababa, Ethiopia

Tel: +251 930 110 004 | +251 937 453 945 | Email: info@anbessaimports.com | anbessa.medical@gmail.com

www.anbessaimports.com

EFDA Registration No. 14280/MDICER/2024