



Anbessa Imports

ANBESSA MEDICAL · STRATEGIC SOURCING & SUPPLIES

40 kW High-Frequency X-Ray System

DIAGNOSTIC IMAGING



40 kW high-frequency digital radiography system with floor-mounted tube column, elevating patient table and wall bucky stand for high-throughput general radiography.

Price: TBC — request a quote

Key Benefits

- Complete radiography room: generator, tube column, table and wall bucky
- Digital flat panel detector for immediate image review
- Built for high daily patient throughput

Specifications

Power	40 kW high-frequency generator
kV Range	40–150 kV
mA Range	10–500 mA
Detector	Flat panel digital (DR)
Configuration	Tube column, table & wall bucky

About Anbessa Imports

Anbessa Imports is the sister company of Anbessa Travel, an Ethiopian business known for its sustainability certification, gold level taxpayer status and 13 years of exceptional customer service. We replicate that same heritage in medical equipment: responsibly sourced from vetted manufacturers, ethically managed through every stage of procurement and compliance, and delivered with installation, training and after-sales support across Ethiopia.

- Sister company of Anbessa Travel
- Sustainability certified, gold level taxpayer
- 13 years of exceptional customer service
- Responsibly sourced and ethically managed

Anbessa Imports — Anbessa Medical

Bole Sub-City, behind the Japanese Embassy, on the road between Bole Sub-City HQ and Egzabherab Church, in front of Chiko Restaurant, Addis Ababa, Ethiopia

Tel: +251 930 110 004 | +251 937 453 945 | Email: info@anbessaimports.com | anbessa.medical@gmail.com

www.anbessaimports.com

EFDA Registration No. 14280/MDICER/2024