



Anbessa Imports

ANBESSA MEDICAL · STRATEGIC SOURCING & SUPPLIES

Biphasic Defibrillator

EMERGENCY



Portable biphasic defibrillator/monitor with AED mode, manual mode, pacing, SpO₂, and 12-lead ECG capability.

Price: TBC — request a quote

Key Benefits

- AED and manual modes with external pacing
- Charges to 200 J in under five seconds
- 200 discharges per battery charge for ambulance use

Specifications

Energy	1–360 J biphasic
Charge Time	< 5 sec to 200J
Display	8.4" color TFT
Battery	200 discharges
Weight	5.8 kg

About Anbessa Imports

Anbessa Imports is the sister company of Anbessa Travel, an Ethiopian business known for its sustainability certification, gold level taxpayer status and 13 years of exceptional customer service. We replicate that same heritage in medical equipment: responsibly sourced from vetted manufacturers, ethically managed through every stage of procurement and compliance, and delivered with installation, training and after-sales support across Ethiopia.

- Sister company of Anbessa Travel
- Sustainability certified, gold level taxpayer
- 13 years of exceptional customer service
- Responsibly sourced and ethically managed

Anbessa Imports — Anbessa Medical

Bole Sub-City, behind the Japanese Embassy, on the road between Bole Sub-City HQ and Egzabherab Church, in front of Chiko Restaurant, Addis Ababa, Ethiopia

Tel: +251 930 110 004 | +251 937 453 945 | Email: info@anbessaimports.com | anbessa.medical@gmail.com

www.anbessaimports.com

EFDA Registration No. 14280/MDICER/2024