



Anbessa Imports

ANBESSA MEDICAL · STRATEGIC SOURCING & SUPPLIES

Genrui KT-42 Hematology Analyzer

LABORATORY

3 PART HEMATOLOGY ANALYZER



3-part differential haematology analyzer using the Coulter principle with RFID closed reagent system, 9 μ L sample volume and wide linearity — ideal for small to medium laboratories.

Price: TBC — request a quote

Key Benefits

- Up to 60 tests per hour with 24 reported parameters
- Only 9 μ L of whole blood per sample
- RFID closed reagent system for consistent results

Specifications

Type	3-part differential
Throughput	Up to 60 tests/hour
Parameters	24 + 3 histograms
Sample Volume	9 μ L whole blood
Reagent System	RFID closed system

About Anbessa Imports

Anbessa Imports is the sister company of Anbessa Travel, an Ethiopian business known for its sustainability certification, gold level taxpayer status and 13 years of exceptional customer service. We replicate that same heritage in medical equipment: responsibly sourced from vetted manufacturers, ethically managed through every stage of procurement and compliance, and delivered with installation, training and after-sales support across Ethiopia.

- Sister company of Anbessa Travel
- Sustainability certified, gold level taxpayer
- 13 years of exceptional customer service
- Responsibly sourced and ethically managed

Anbessa Imports — Anbessa Medical

Bole Sub-City, behind the Japanese Embassy, on the road between Bole Sub-City HQ and Egzabherab Church, in front of Chiko Restaurant, Addis Ababa, Ethiopia

Tel: +251 930 110 004 | +251 937 453 945 | Email: info@anbessaimports.com | anbessa.medical@gmail.com

www.anbessaimports.com

EFDA Registration No. 14280/MDICER/2024