



Anbessa Imports

ANBESSA MEDICAL · STRATEGIC SOURCING & SUPPLIES

Mobile C-Arm Fluoroscopy System

DIAGNOSTIC IMAGING



Mobile C-arm imaging system with high-resolution fluoroscopy, dual-monitor workstation and motorised positioning for surgical, orthopaedic and interventional procedures.

Price: TBC — request a quote

Key Benefits

- Intra-operative imaging for orthopaedic and interventional theatres
- Motorised positioning with $\pm 180^\circ$ orbital rotation
- Dual-monitor workstation for surgeon and radiographer

Specifications

Generator	Up to 5 kW
Imaging	Fluoroscopy, pulsed & digital radiography
Detector	9"/12" image intensifier or flat panel
Monitors	Dual high-resolution workstation
Rotation	$\pm 180^\circ$ orbital, motorised lift

About Anbessa Imports

Anbessa Imports is the sister company of Anbessa Travel, an Ethiopian business known for its sustainability certification, gold level taxpayer status and 13 years of exceptional customer service. We replicate that same heritage in medical equipment: responsibly sourced from vetted manufacturers, ethically managed through every stage of procurement and compliance, and delivered with installation, training and after-sales support across Ethiopia.

- Sister company of Anbessa Travel
- Sustainability certified, gold level taxpayer
- 13 years of exceptional customer service
- Responsibly sourced and ethically managed

Anbessa Imports — Anbessa Medical

Bole Sub-City, behind the Japanese Embassy, on the road between Bole Sub-City HQ and Egzabherab Church, in front of Chiko Restaurant, Addis Ababa, Ethiopia

Tel: +251 930 110 004 | +251 937 453 945 | Email: info@anbessaimports.com | anbessa.medical@gmail.com

www.anbessaimports.com

EFDA Registration No. 14280/MDICER/2024