



Anbessa Imports

ANBESSA MEDICAL · STRATEGIC SOURCING & SUPPLIES

Surgical Operating Table (Pediatric)

SURGICAL EQUIPMENT



Specialized pediatric surgical table with precise height and tilt adjustments, radiolucent tabletop, and child-safe side rails for safe pediatric procedures.

Price: TBC — request a quote

Key Benefits

- Radiolucent carbon-fibre tabletop for intra-operative imaging
- Precise height, lateral tilt and Trendelenburg adjustment
- Paediatric restraints and child-safe side rails included

Specifications

Load Capacity	120 kg
Height Range	600–1000 mm
Tabletop	Radiolucent carbon fiber
Tilt	$\pm 25^{\circ}$ lateral, $\pm 30^{\circ}$ Trendelenburg
Accessories	Pediatric restraints included

About Anbessa Imports

Anbessa Imports is the sister company of Anbessa Travel, an Ethiopian business known for its sustainability certification, gold level taxpayer status and 13 years of exceptional customer service. We replicate that same heritage in medical equipment: responsibly sourced from vetted manufacturers, ethically managed through every stage of procurement and compliance, and delivered with installation, training and after-sales support across Ethiopia.

- Sister company of Anbessa Travel
- Sustainability certified, gold level taxpayer
- 13 years of exceptional customer service
- Responsibly sourced and ethically managed

Anbessa Imports — Anbessa Medical

Bole Sub-City, behind the Japanese Embassy, on the road between Bole Sub-City HQ and Egzabherab Church, in front of Chiko Restaurant, Addis Ababa, Ethiopia

Tel: +251 930 110 004 | +251 937 453 945 | Email: info@anbessaimports.com | anbessa.medical@gmail.com

www.anbessaimports.com

EFDA Registration No. 14280/MDICER/2024